



Oxivir Excel Foam Mini Guide

Available in 750ml RTU foaming trigger



Use for disinfecting:



Non invasive utensils
and equipment



Washable surfaces

Do not use on marble



Ready to use



Fully bactericidal and Virucidal in
30 seconds including gram-negative
bacteria and non-enveloped viruses




Using...

1



Remove any gross
soilage and debris
with pre-soaked
paper towels/
disposable cloths.

2



Spray product onto
clean surface.

3



Allow 30 second
contact time.

4



Wipe away
excess.

5



If dealing with
a body spillage,
dispose as
contaminated
waste.

