



TASKI® Sani Clonet

Toilet cleaner

Description

Ready to use, heavy-duty cleaner for the effective removal of limescale and other deposits from toilet bowls and urinals.

Key properties

- Effective combination of three different acids
- Clings to sloped surfaces for longer contact time
- Has an angled neck bottle to reach under the toilet rim easily
- Unique, patented Odour Neutralizing Technology (O.N.T.)

Benefits

- Quickly dissolves calcium/lime deposits
- Long-lasting action on vertical/sloped surfaces
- Easily reaches hidden areas e.g. under the toilet rim

Use instructions

Dosage:

Ready to use, no dilution required.

Application:

1. Flush the toilet/urinal.
2. Apply the product under the rim, around the bowl, onto water/flush paths and water outlets.
3. Allow the product to act for 5-10 minutes.
4. Scrub with a brush and flush to rinse.
5. Increase acting time (1-2 hours) in case of heavy calcium/lime build-up.



© A.I.S.E.



TASKI® Sani Clonet

Technical data

Appearance: Clear, red liquid

pH value (neat): 0.7

pH value (in use): 0.7

Relative density (20°C): 1.06 g/cm³

The above data is typical of normal production and should not be taken as a specification.

Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet; sds.diversey.com. Store in original closed containers away from extremes of temperature. Store in original closed containers away from extremes of temperature..

Product compatibility

Do not mix the product with chlorinated products and do not use on acid sensitive surfaces (e.g. marble, enamel, grouting). Immediately remove any drops or splashes from acid-sensitive surfaces and rinse thoroughly with water. Can be used on intact porcelain and glazed sanitary installations. Test material compatibility in a small, inconspicuous place before use.

Environmental information

The surfactants used in this product are biodegradable in compliance with the requirement of EU Detergent Regulation, EC 648/2004.

Available pack sizes

TASKI Sani Clonet is available in 6x0.75L, 2x5L and 1x10L